

# NeO<sub>2</sub>-Safe™

## Provides access without disconnecting circuits



NeO<sub>2</sub>-Safe Infant Ventilator Circuit Manifold provides access to the endotracheal tube without disrupting ventilation, PEEP or CPAP. NeO<sub>2</sub>-Safe fits all neonatal circuits and is manufactured from latex free, hypoallergenic materials. The NeO<sub>2</sub>-Safe access port maintains a closed airway circuit system by sealing the valve around the catheter. The access port is a Duckbill valve made of medical grade silicone.

When infants are disconnected from the ventilator for specific procedures, i.e., suctioning, many infants rapidly desaturate due to the loss of positive end-expiratory pressure (PEEP) and their inability to maintain FRC. A method to allow for suctioning without disconnect from the ventilator circuit may prevent or minimize this occurrence. The NeO<sub>2</sub>-Safe used in-line with a 10 mm ventilator circuit allows access for suctioning without disconnect from the ventilator, effectively reducing the potential for hypoxemia associated with open circuit suctioning. The NeO<sub>2</sub>-Safe maintains PEEP or CPAP during suction procedures.

The NeO<sub>2</sub>-Safe is intended for single patient use application and should not be reused.

### Ordering Information

11100	NeO <sub>2</sub> -Safe (100/case)
11103	NeO <sub>2</sub> -Safe Cap (100/case)
11105	NeO <sub>2</sub> -Safe (50/box)

### Complementary Products

11700	Hybrid Baby Tape™ (100/case)
11800	Hybrid Pediatric Tape™ (100/case)

### B&B Products are available from

B&B Medical Technologies and finer specialty care distributors.

Visit [www.BandB-Medical.com](http://www.BandB-Medical.com) or contact us today at 1.800.242.8778.

### SAFE

Enhanced durability of the Duckbill valve to maintain ventilation and PEEP

Latex free

Hypoallergenic

### COST EFFECTIVE

Dual purpose manifold (suction and ventilator circuit interface)

Reduction of costs that can be associated with closed suction catheter systems

### CONVENIENT

Adapts to infant ventilator circuits (10 mm/reusable or single patient use)

Easy access to the airway for suction procedure via the Duckbill valve

### VERSATILE

Designed for use in the NICU, Labor and Delivery and transport situations

11100C2

© 2009 B&B Medical Technologies. All rights reserved.

"A World of Products for Better Breathing" is a service mark of B&B Medical Technologies.

Toll-free: +1.800.242.8778

Tel: +1.760.929.9972

Fax: +1.760.929.9953

2734 Loker Avenue West, Suite M  
Carlsbad, CA 92010 USA

