

OPTI-MIST II/OPTI-MIST CLEAR Nebulizers & Kits

Improved OPTI-MIST II and New OPTI-MIST CLEAR are two types of small volume nebulizers with easy to use designs. They are spill proof and provide fast and effective medication delivery. Medication capacity is 5cc/ml and is conveniently graduated on medication cups. The OPTI-MIST style locks with a simple 1/4 turn.



Ref#750-E



Ref#753-E



Ref#774-E



Ref#752-E

NEBULIZER SPECIFICATION CHART

Criteria / Product	Requirement	OPTI-MIST	OPTI-MIST CLEAR	UP-MIST
TEST AIRFLOW		6 LPM	6LPM	6 LPM
MMAD	1.0-5.0 micron	3.00	2.60	3.60
Minimum Aerosol Output (cc/ml) after 10 Min. Operation	Min. 1.500 cc/ml	1.90	2.60	2.10
Minimum Fill Volume (cc/ml) for 10 Min. Operation	Lowest Possible	2.00	1.60	2.30
Residual Volume (cc/ml) 5 seconds after mist ends	Lowest Possible	1.10	1.50	2.60

All results tested from vertical position

ALL PRODUCTS LATEX-FREE

OPTI-MIST II NEBULIZER & KITS

* Ref#	Description	Qty/Cs
750-E	Nebulizer Only	NEW 50
3 752-E	Nebulizer, 7 ft. (2.1 m) Oxygen Tubing, "T" Piece, Mouth Piece	NEW 50
3 753-E	Nebulizer, 7 ft. (2.1 m) Oxygen Tubing, Adult Aerosol Mask, Elongated	NEW 50
758-E	Nebulizer, No Tubing, "T" Piece, Mouth Piece	NEW 50
3 759-E	Nebulizer, 7 ft. (2.1 m) Oxygen Tubing, "T" Piece, Mouth Piece, 6 in. (15 cm) Aerosol Tubing	NEW 50
3 760-E	Nebulizer, 7 ft. (2.1 m) Oxygen Tubing, "T" Piece (22mm O.D./22mm I.D.), Mouth Piece	NEW 50
3 229-E	Nebulizer, 7 ft. (2.1 m) Oxygen Tubing, Pediatric Aerosol Mask, Elongated	NEW 50

OPTI-MIST CLEAR NEBULIZER & KITS

* Ref#	Description	Qty/Cs
774-E	Nebulizer Only	NEW 50
3 775-E	Nebulizer, 7 ft. (2.1 m) Oxygen Tubing, "T" Piece, Mouth Piece, 6 in. (15 cm) Aerosol Tubing	NEW 50
3 776-E	Nebulizer, 7 ft. (2.1 m) Oxygen Tubing, Adult Aerosol Mask, Elongated	NEW 50
3 777-E	Nebulizer, 7 ft. (2.1 m) Oxygen Tubing, Pediatric Aerosol Mask, Elongated	NEW 50
778-E	Nebulizer, 7 ft. (2.1 m) Oxygen Tubing, "T" Piece, Mouth Piece	NEW 50

*For FITS-ALL Oxygen Connector add the digit "3" in front of Ref# (e.g. 3752-E)

Source:Nelson Lab