

Hyperinflation Sets

TOTAL CONTROL AT YOUR FINGERTIPS

Complete
Control



Westmed

When Your Care is Critical

Hyperinflation Sets

TOTAL CONTROL AT YOUR FINGERTIPS

Westmed hyperinflation sets are precision engineered for versatility and functionality. They are the perfect choice for the NICU or post-operative transport.

- The high profile, knurled adjustment knob places full pressure control at the clinician's fingertips.
- ISO dimensions to ensure a secure, leak resistant attachment.
- Precision molded bleed orifice minimizes inadvertent PEEP or CPAP, while providing safe, effective pressure-limitation.
- Configurations available with disposable manometers as depicted in the photo above.
- Default safety bleed aids in prevention of over inflation.



Transport Sets		Case Quantity
0413	Transport Elbow, 8' Oxygen Tubing w/Threaded Nut	20
6214	Norman Elbow, 7' Oxygen Tubing w/Threaded Nut, Cut-Off Tail	20
8073	Transport Circuit, 7' Oxygen Tubing w/Threaded Nut	20
7-N7910	Hyperinflation Elbow, 7' Oxygen Tubing w/Threaded Nut, 7' Pressure Line	50
Hyperinflation Sets w/Oxygen and Pressure Tubing		
7-N7905-62	½ Liter Latex Free Bag, 84" O2 & Pressure Lines	12
7-N7901-62	1 Liter Latex Free Bag, 84" O2 & Pressure Lines	12
7-N7902-62	2 Liter Latex Free Bag, 84" O2 & Pressure Lines	12
7-N7903-62	3 Liter Latex Free Bag, 84" O2 & Pressure Lines	12
Hyperinflation Sets w/Masks		
7-N7905-I	½ Liter Latex Free Bag, 84" O2 & Pressure Lines, Infant Mask	12
7-N7905-N	½ Liter Latex Free Bag, 84" O2 & Pressure Lines, Neonatal Mask	12
Hyperinflation Sets w/Disposable Manometers		
7-N7905-MAN	½ Liter Latex Free Bag, 84" O2 Line, Manometer	12
7-N7901-MAN	1 Liter Latex Free Bag, 84" O2 Line, Manometer	12
Hyperinflation Sets w/Disposable Manometers and Masks		
7-N7905-I-MAN	½ Liter Latex Free Bag, 84" O2 Line, Infant Mask, Manometer	12
7-N7905-N-MAN	½ Liter Latex Free Bag, 84" O2 Line, Neonatal Mask, Manometer	12
Hyperinflation Accessories		
0497	Manometer Cover	1
0499	Manometer	1



Please log onto www.westmedinc.com for a complete listing of Westmed anesthesia and respiratory products.

© 2007 Westmed, Inc. All rights reserved. PN 71408-09

Westmed

When Your Care is Critical

Westmed, Inc.

5580 S. Nogales Hwy.,
Tucson, AZ 85706

Phone: 800-975-7987

Fax: 520-294-6061

www.westmedinc.com