


KING VALUE MASK



The King Value Mask is clinically proven to provide an efficient seal during anesthesia. The King Value Mask has a thin wall cushion that conforms to the patient’s anatomy. The mask crown is easy to hold and assists in creating a good seal. The clear crown and cushion allows the clinician to easily see the patient’s breath. Sizes 3, 4, 5, and 6 come with a hook ring to allow easy attachment of a head strap.

Available in six sizes from neonate to large adult.

Size	Patient	Product Code	Quantity
Size 1	Neonate	1915	50 per case
Size 2	Infant	1925	50 per case
Size 3	Toddler	1935	50 per case
Size 4	Small Adult	1945	50 per case
Size 5	Medium Adult	1955	50 per case
Size 6	Larger Adult	1965	50 per case