



Ambu aScope™ 4 RhinoLaryngo

**Guaranteed sterility**  
**Single-use efficiency**



# Single-use – the simple way to improve productivity and safety

## Guaranteed sterility

When it comes to disinfection techniques for reusable endoscopes, guidelines within ENT are becoming increasingly strict. As facilities contend with these more rigorous cleaning requirements, there is a growing need for a better way to ensure safety without compromising efficiency.

---

Introducing aScope 4  
RhinoLaryngo, the single-use  
flexible video endoscope that  
is always sterile with no risk  
of cross-contamination.



# A hassle-free solution for a simplified workflow

## Single-use efficiency

aScope 4 RhinoLaryngo allows you to simplify your workflow and increase productivity – for the benefit of both patients and staff.

Sterile straight from the pack, aScope 4 RhinoLaryngo eliminates a number of costly and time-consuming tasks associated with reusable scopes, including cleaning and repairs. As it is single-use, there is no need to establish a new, more complicated workflow or make significant investments due to increased cleaning requirements.

## Single-use transparency

aScope 4 RhinoLaryngo also offers the assurance of complete cost transparency thanks to a simple per unit cost, which also helps you avoid the extra expenses and delays that come with repairs and reprocessing.



---

# High quality imaging and intuitive design

## Clear and sharp images

aScope 4 RhinoLaryngo gives you clear visibility of the upper airway, and the sharp images and adaptive light control enable easy identification of anatomical structures. aScope 4 RhinoLaryngo also features high bending angles and precise tip motion to make navigation in the upper airway easy.

## Lightweight comfort

No matter how many procedures you do each day, aScope 4 RhinoLaryngo is always comfortable to hold thanks to its lightweight design. The handle has a rippled surface and an ergonomic grip, and the design makes it intuitive to use and easy to manoeuvre throughout each examination.

With its external monitor, aScope 4 RhinoLaryngo also promotes a better working posture compared to traditional fibrescopes with eyepieces, and it allows for a more pleasant examination experience for patients.




# Portable video imaging with aView™

aScope 4 RhinoLaryngo connects to aView, our high-resolution monitor that is portable, simple to set up and intuitive to use. By providing high quality images on an external monitor, aView supports training, knowledge sharing and patient involvement. aView also enables you to record and store images and videos of procedures for later reference.

The compact aView monitor fits easily into small office spaces. It can be quickly mounted on an IV pole or transported by hand for additional convenience.

As an integrated part of our single-use portfolio, aView is compatible with all single-use flexible scopes from Ambu – including aScope 4 RhinoLaryngo Slim and Intervention.



# A complete endoscope solution for ENT

## Hassle-free safety

Working channels present particular problems when it comes to cleaning reusable scopes. These scopes require full reprocessing, and effective disinfection is complicated by their long, narrow lumens and delicate materials. Even the most stringent cleaning procedures do not guarantee the complete removal of biofilm from a reusable scope.

Sterile straight from the pack, single-use aScope 4 RhinoLaryngo Slim and Intervention ensure patient safety by eliminating the risk of device-related cross-contamination.

The aScope 4 RhinoLaryngo portfolio includes two single-use ENT endoscopes to cover your needs in flexible nasopharyngoscopy procedures.

- aScope 4 RhinoLaryngo Slim\* has a small outer diameter of 3.0 mm to help minimize patient discomfort during procedures, and its precise tip motion and 130-degree bending angle support easy manoeuvrability and navigation in the upper airway. The Slim scope has a length of 300 mm.
- With its 2.2 mm working channel, the aScope 4 RhinoLaryngo Intervention\*\* is designed for examinations such as biopsies, foreign body removal and suctioning procedures. The Intervention scope offers the same manoeuvrability as the Slim scope, and has a length of 350 mm and an outer diameter of 5.0 mm.



**Ambu Inc.**

6230 Old Dobbin Lane  
Suite 250  
Columbia, MD 21045  
800 262 8462

For more information  
**Visit [ambuUSA.com](http://ambuUSA.com)**

Distributed By:  **CAREstream Medical**

Unit 103-19099 25 Avenue, Surrey, BC V3Z 3V2

Tel: 1.888.310.2186 Fax: 1.888.310.2187

[info@carestreammedical.com](mailto:info@carestreammedical.com) [www.carestreammedical.com](http://www.carestreammedical.com)

**Trusted innovation from Ambu**

Since 1937 Ambu has been bringing breakthrough healthcare solutions to life within anaesthesia, patient monitoring and diagnostics, and emergency care. Today, we are dedicated to improving safety for millions of patients and healthcare professionals around the world through the advancement of single-use devices, including aScope – the world's first single-use flexible video scope.

Our commitment to bringing new ideas, exceptional product performance and superior service to our customers has made Ambu one of the most recognised medical companies in the world. Headquartered near Copenhagen, Denmark, Ambu employs around 2,500 people across Europe, North America and Asia Pacific. For more information, please visit [www.ambu.com](http://www.ambu.com)

\* aScope 4 RhinoLaryngo Slim is pending registrations and is not released for sale.

\*\* aScope 4 RhinoLaryngo Intervention is not available for sale in all countries.  
Please contact your local Ambu sales representative.

US: Rx only