

# Canada Product List

Effective March 4, 2013

# RESMED

## Lasting Results

QUALITY SYSTEMS

COMPLETE SOLUTIONS

RELIABLE PRODUCTS

• Swift FX and Swift FX Bella Gray Nasal Pillows Systems	2
• Swift FX for Her and Swift FX Bella Nasal Pillows Systems	3
• Swift LT and Swift LT for Her Nasal Pillows Systems	4,5
• Mirage Swift II Nasal Pillows System and Accessories	6
• Mirage FX and Mirage FX for Her Nasal Mask Systems	7
• Mirage SoftGel Nasal Mask System and Accessories	8,9
• Mirage Activa LT Nasal Mask System and Accessories	10,11
• Mirage SoftGel/Mirage Activa LT ConvertAble Pack	12
• Mirage Activa Nasal Mask Systems and Accessories	13
• Mirage Micro and Mirage Micro for Kids Nasal Mask Systems	14,15
• Ultra Mirage II Nasal Mask System and Accessories	16
• Pixi Pediatric Nasal Mask System	17
• Mirage Vista and Mirage Kidsta Nasal Mask Systems	18
• Infant Nasal Mask System and Accessories	19
• Quattro FX and Quattro FX for Her Full Face Mask Systems	20,21
• Mirage Quattro Full Face Mask System and Accessories	22,23
• Ultra Mirage Full Face Mask System and Accessories	24,25
• Mirage Liberty Full Face Mask System and Accessories	26,27
• Mask Spares and Accessories	28
• S9 APAP, S9 CPAP and S9 VPAP Flow Generators	29,30,31
• S9 Humidifiers and Accessories	32
• Respiratory Care Products	33,34
• Software, Modules and Accessories	35,36
• ApneaLink/ApneaLink Plus and Accessories	37,38
• Flow Generator Accessories and Filters	39,40,41
• Heated Humidification Systems and Accessories	42



Item #	Description
61500	<b>Swift™ FX nasal pillows system</b> Includes small, medium, and large size dual-wall pillows, and headgear
61568	<b>Swift™ FX Bella Gray nasal pillows system</b> Includes small, medium and large size dual-wall pillows, Bella Gray Loops accessory and Swift FX headgear with soft wraps
PILLOWS	
61520	<b>Extra small dual-wall pillow</b>
61521	<b>Small dual-wall pillow</b>
61522	<b>Medium dual-wall pillow</b>
61523	<b>Large dual-wall pillow</b>
FRAME	
61510	<b>Swift FX frame system with extra small dual-wall pillows - no headgear</b>
61511	<b>Swift FX frame system with small dual-wall pillows - no headgear</b>
61512	<b>Swift FX frame system with medium dual-wall pillows - no headgear</b>
61513	<b>Swift FX frame system with large dual-wall pillows - no headgear</b>
61528	<b>Swift FX frame assembly - no dual-wall pillows, no headgear</b>
HEADGEAR	
61529	<b>Swift FX headgear assembly</b>
61583	<b>Swift FX headgear assembly with Bella Gray Loops accessory (headgear, gray and light blue back straps, gray soft wraps, Bella Gray Loops)</b>
ACCESSORIES	
61530	<b>Soft wraps - 2 pack</b>
61121	<b>Distal pressure port attachment</b>
60549	<b>Pressure port plug - 1 pack</b>
60550	<b>Pressure port plug - 10 pack</b>
61582	<b>Swift FX Bella Gray Loops – 2 pack (gray)</b>

**Swift™ FX for Her  
Swift™ FX Bella  
Nasal Pillows System**



Item #	Description
61540	<b>Swift™ FX for Her nasal pillows system</b> Includes extra small, small, and medium size dual-wall pillows, and headgear
61560	<b>Swift™ FX Bella nasal pillows system</b> Includes extra small, small, and medium size dual-wall pillows, Bella Loops accessory and Swift FX for Her headgear with soft wraps
<b>PILLOWS</b>	
61520	<b>Extra small dual-wall pillow</b>
61521	<b>Small dual-wall pillow</b>
61522	<b>Medium dual-wall pillow</b>
61523	<b>Large dual-wall pillow</b>
<b>FRAME</b>	
61510	<b>Swift FX frame system with extra small dual-wall pillows - no headgear</b>
61511	<b>Swift FX frame system with small dual-wall pillows - no headgear</b>
61512	<b>Swift FX frame system with medium dual-wall pillows - no headgear</b>
61513	<b>Swift FX frame system with large dual-wall pillows - no headgear</b>
61528	<b>Swift FX frame assembly - no dual-wall pillows, no headgear</b>
<b>HEADGEAR</b>	
61543	<b>Swift FX for Her headgear assembly (headgear, pink and light blue back straps, pink and gray soft wraps)</b>
61580	<b>Swift FX for Her headgear assembly with Bella Loops accessory (headgear, pink and light blue back straps, pink and gray soft wraps, Bella Loops)</b>
<b>ACCESSORIES</b>	
61530	<b>Soft wraps - 2 pack (gray)</b>
61544	<b>Soft wraps - 2 pack (pink)</b>
61121	<b>Distal pressure port attachment</b>
60549	<b>Pressure port plug - 1 pack</b>
60550	<b>Pressure port plug - 10 pack</b>
61581	<b>Swift FX Bella Loops – 2 pack (pink)</b>



Item #	Description
60560	<b>Swift™ LT nasal pillows system</b> Includes small, medium, and large size dual-wall pillows, and headgear
<b>PILLOWS</b>	
60574	<b>Extra small dual-wall pillow</b>
60571	<b>Small dual-wall pillow</b>
60572	<b>Medium dual-wall pillow</b>
60573	<b>Large dual-wall pillow</b>
<b>FRAME</b>	
60585	<b>Swift LT frame system with extra small dual-wall pillows - no headgear</b>
60582	<b>Swift LT frame system with small dual-wall pillows - no headgear</b>
60583	<b>Swift LT frame system with medium dual-wall pillows - no headgear</b>
60584	<b>Swift LT frame system with large dual-wall pillows - no headgear</b>
60581	<b>Swift LT frame assembly - no dual-wall pillows, no headgear</b>
60575	<b>Frame only</b>
<b>HEADGEAR</b>	
60578	<b>Headgear assembly</b>
<b>ACCESSORIES</b>	
61121	<b>Distal pressure port attachment</b>
60558	<b>Headgear buckle (back) - 1 pack</b>
60559	<b>Headgear buckle (back) - 10 pack</b>
60579	<b>Headgear buckle (top) - 1 pack</b>
60580	<b>Headgear buckle (top) - 10 pack</b>
60577	<b>Short tube assembly, includes elbow, short tube and swivel - 1 pack</b>
60742	<b>Short tube assembly, includes elbow, short tube and swivel - 10 pack</b>
60586	<b>Tube retainer (plastic)</b>
60557	<b>Tube retainer (velcro)</b>
60549	<b>Pressure port plug - 1 pack</b>
60550	<b>Pressure port plug - 10 pack</b>
60599	<b>Soft wraps - 2 pack</b>



Item #	Description
60588	<b>Swift™ LT for Her nasal pillows system</b> Includes extra small, small and medium size dual-wall pillows, and headgear
<b>PILLOWS</b>	
60574	<b>Extra small dual-wall pillow</b>
60571	<b>Small dual-wall pillow</b>
60572	<b>Medium dual-wall pillow</b>
60573	<b>Large dual-wall pillow</b>
<b>FRAME</b>	
60585	<b>Swift LT frame system with extra small dual-wall pillows - no headgear</b>
60582	<b>Swift LT frame system with small dual-wall pillows - no headgear</b>
60583	<b>Swift LT frame system with medium dual-wall pillows - no headgear</b>
60584	<b>Swift LT frame system with large dual-wall pillows - no headgear</b>
60581	<b>Swift LT frame assembly - no dual-wall pillows, no headgear</b>
60575	<b>Frame only</b>
<b>HEADGEAR</b>	
60595	<b>Headgear assembly - light blue</b>
<b>ACCESSORIES</b>	
61121	<b>Distal pressure port attachment</b>
60740	<b>Headgear buckle (back) - 1 pack - light blue</b>
60741	<b>Headgear buckle (back) - 10 pack - light blue</b>
60597	<b>Headgear buckle (top) - 1 pack - light blue</b>
60596	<b>Headgear buckle (top) - 10 pack - light blue</b>
60577	<b>Short tube assembly, includes elbow, short tube and swivel - 1 pack</b>
60742	<b>Short tube assembly, includes elbow, short tube and swivel - 10 pack</b>
60599	<b>Soft wraps - 2 pack</b>
60598	<b>Tube retainer (plastic) - light blue</b>
60549	<b>Pressure port plug - 1 pack</b>
60550	<b>Pressure port plug - 10 pack</b>



Item #	Description
60512	<b>Mirage Swift™ II nasal pillows system</b> Includes small, medium, and large size pillow sleeves, and headgear
<b>PILLOW SLEEVES</b>	
60541	<b>Small pillow sleeve</b>
60542	<b>Medium pillow sleeve</b>
60543	<b>Large pillow sleeve</b>
<b>FRAME</b>	
60537	<b>Mirage Swift II frame system with small pillows - no headgear</b>
60538	<b>Mirage Swift II frame system with medium pillows - no headgear</b>
60539	<b>Mirage Swift II frame system with large pillows - no headgear</b>
60544	<b>Mirage Swift II frame assembly - no pillows, no headgear</b>
60536	<b>Frame with clip and cap only</b>
<b>HEADGEAR</b>	
60526	<b>Headgear</b>
<b>ACCESSORIES</b>	
61121	<b>Distal pressure port attachment</b>
60523	<b>Headgear buckle – 4 pack</b>
60540	<b>Frame clip</b>
60524	<b>Seal ring - 10 pack</b>
60525	<b>Seal cap</b>
60528	<b>Short tube assembly, includes elbow, short tube and swivel - 1 pack</b>
60529	<b>Short tube assembly, includes elbow, short tube and swivel - 10 pack</b>
60530	<b>Tube retainer - 3 pack</b>
60549	<b>Pressure port plug - 1 pack</b>
60550	<b>Pressure port plug - 10 pack</b>



Item #	Description
62103	<b>Mirage™ FX nasal mask complete system</b> Includes cushion and headgear
62118	<b>Mirage FX nasal mask complete system - wide</b> Includes cushion and headgear
62128	<b>Mirage™ FX for Her nasal mask complete system</b> Includes cushion and headgear
62109	<b>Mirage FX for Her nasal mask complete system - small</b> Includes small cushion and small headgear
<b>CUSHIONS</b>	
62111	<b>Mirage FX and Mirage FX for Her cushion - standard</b>
62125	<b>Mirage FX cushion - wide</b>
62136	<b>Mirage FX cushion - small</b>
<b>FRAME SYSTEM</b>	
62112	<b>Mirage FX and Mirage FX for Her frame system with cushion - standard</b>
62126	<b>Mirage FX frame system with cushion - wide</b>
62137	<b>Mirage FX frame system with cushion - small</b>
<b>HEADGEAR</b>	
62138	<b>Mirage FX headgear - small (gray)</b>
62110	<b>Mirage FX headgear</b>
62129	<b>Mirage FX for Her headgear</b>



Item #	Description
61600	<b>Mirage™ SoftGel nasal mask complete system - small</b> Includes cushion and headgear
61601	<b>Mirage SoftGel nasal mask complete system - medium</b> Includes cushion and headgear
61602	<b>Mirage SoftGel nasal mask complete system - large</b> Includes cushion and headgear
61603	<b>Mirage SoftGel nasal mask complete system - large wide</b> Includes cushion and headgear
<b>CUSHIONS</b>	
61631	<b>Small cushion and clip</b>
61632	<b>Medium cushion and clip</b>
61633	<b>Large cushion and clip</b>
61634	<b>Large wide cushion and clip</b>
<b>FRAME</b>	
61626	<b>Mirage SoftGel nasal mask frame system with small cushion - no headgear, includes headgear clips</b>
61627	<b>Mirage SoftGel nasal mask frame system with medium cushion - no headgear, includes headgear clips</b>
61628	<b>Mirage SoftGel nasal mask frame system with large cushion - no headgear, includes headgear clips</b>
61629	<b>Mirage SoftGel nasal mask frame system with large wide cushion - no headgear, includes headgear clips</b>
60196	<b>Mirage Activa LT/SoftGel nasal mask complete frame assembly - small, no cushion, no headgear, includes headgear clips</b>
60175	<b>Mirage Activa LT/SoftGel nasal mask complete frame assembly - standard, no cushion, no headgear, includes headgear clips</b>
60180	<b>Mask frame only - standard</b>





<b>HEADGEAR</b>	
16117	Universal headgear for Ultra Mirage II, Micro, Activa, Activa LT, SoftGel, and Quattro - medium
16118	Universal headgear for Ultra Mirage II, Micro, Micro for Kids, Activa, Activa LT, SoftGel, Ultra Mirage Full Face Mask and Quattro - small
16119	Universal headgear for Ultra Mirage II, Micro, Activa, Activa LT, SoftGel, Ultra Mirage Full Face Mask and Quattro- large
<b>ACCESSORIES</b>	
16399	Elbow assembly - includes elbow and swivel
60176	Cushion clip - standard
60197	Cushion clip - small
60123	Forehead pad - 1 pack
60124	Forehead pad - 10 pack
16393	Forehead support
16569	Headgear clip - 2 pack
16734	Headgear clip - 10 pack
16570	Port cap - 2 pack
16565	Quick release swivel
61289	Forehead Dial
16387	Elbow



Item #	Description
60182	<b>Mirage Activa™ LT nasal mask complete system - small</b> Includes cushion and headgear
60148	<b>Mirage Activa LT nasal mask complete system - medium</b> Includes cushion and headgear
60149	<b>Mirage Activa LT nasal mask complete system - large</b> Includes cushion and headgear
60150	<b>Mirage Activa LT nasal mask complete system - large wide</b> Includes cushion and headgear
<b>CUSHIONS</b>	
60198	<b>Small cushion and clip</b>
60177	<b>Medium cushion and clip</b>
60178	<b>Large cushion and clip</b>
60179	<b>Large wide cushion and clip</b>
<b>FRAME</b>	
60195	<b>Mirage Activa LT nasal mask frame system with small cushion - no headgear, includes headgear clips</b>
60172	<b>Mirage Activa LT nasal mask frame system with medium cushion - no headgear, includes headgear clips</b>
60173	<b>Mirage Activa LT nasal mask frame system with large cushion - no headgear, includes headgear clips</b>
60174	<b>Mirage Activa LT nasal mask frame system with large wide cushion - no headgear, includes headgear clips</b>
60196	<b>Mirage Activa LT/SoftGel nasal mask complete frame assembly - small, no cushion, no headgear, includes headgear clips</b>
60175	<b>Mirage Activa LT/SoftGel nasal mask complete frame assembly - standard, no cushion, no headgear, includes headgear clips</b>
60193	<b>Mask frame only - small</b>
60180	<b>Mask frame only - standard</b>



<b>HEADGEAR</b>	
16117	Universal headgear for Ultra Mirage II, Micro, Activa, Activa LT, SoftGel, and Quattro - medium
16118	Universal headgear for Ultra Mirage II, Micro, Micro for Kids, Activa, Activa LT, SoftGel, Ultra Mirage Full Face Mask and Quattro - small
16119	Universal headgear for Ultra Mirage II, Micro, Activa, Activa LT, SoftGel, Ultra Mirage Full Face Mask and Quattro- large
<b>ACCESSORIES</b>	
16399	Elbow assembly - includes elbow and swivel
60176	Cushion clip - standard
60197	Cushion clip - small
60123	Forehead pad - 1 pack
60124	Forehead pad - 10 pack
16393	Forehead support
16569	Headgear clip - 2 pack
16734	Headgear clip - 10 pack
16570	Port cap - 2 pack
16565	Quick release swivel
61289	Forehead Dial
16387	Elbow

Item #	Description
61604	<b>Mirage SoftGel/Mirage Activa LT ConvertAble Pack - small</b> Includes one SoftGel mask, one Activa LT cushion and headgear
61609	<b>Mirage SoftGel/Mirage Activa LT ConvertAble Pack - medium</b> Includes one SoftGel mask, one Activa LT cushion and headgear
61615	<b>Mirage SoftGel/Mirage Activa LT ConvertAble Pack - large</b> Includes one SoftGel mask, one Activa LT cushion and headgear
61620	<b>Mirage SoftGel/Mirage Activa LT ConvertAble Pack - large wide</b> Includes one SoftGel mask, one Activa LT cushion and headgear



Item #	Description
60100	<b>Mirage Activa™ nasal mask complete system - standard</b> Includes cushion and headgear
60101	<b>Mirage Activa™ nasal mask complete system - large</b> Includes cushion and headgear
60102	<b>Mirage Activa™ nasal mask complete system - shallow</b> Includes cushion and headgear
<b>CUSHIONS</b>	
60117	<b>Standard cushion</b>
60118	<b>Large cushion</b>
60119	<b>Shallow cushion</b>
<b>FRAME</b>	
60111	<b>Mirage Activa nasal mask frame system with standard cushion - no headgear, includes headgear clips</b>
60112	<b>Mirage Activa nasal mask frame system with large cushion - no headgear, includes headgear clips</b>
60113	<b>Mirage Activa nasal mask frame system with shallow cushion - no headgear, includes headgear clips</b>
60144	<b>Mirage Activa nasal mask complete frame assembly - no cushion, no headgear, includes headgear clips</b>
60121	<b>Mask frame - one size fits all, mask frame only</b>
<b>HEADGEAR</b>	
16117	<b>Universal headgear for Ultra Mirage II, Micro, Activa, Activa LT, SoftGel, and Quattro - medium</b>
16118	<b>Universal headgear for Ultra Mirage II, Micro, Micro for Kids, Activa, Activa LT, SoftGel, Ultra Mirage Full Face Mask and Quattro - small</b>
16119	<b>Universal headgear for Ultra Mirage II, Micro, Activa, Activa LT, SoftGel, Ultra Mirage Full Face Mask and Quattro - large</b>
<b>ACCESSORIES</b>	
60120	<b>Cushion clip</b>
60125	<b>Elbow assembly - includes exhaust vent</b>
60123	<b>Forehead pad - 1 pack</b>
60124	<b>Forehead pad - 10 pack</b>
60122	<b>Forehead support</b>
60115	<b>Headgear clips - 2 pack</b>
60116	<b>Headgear clips - 10 pack</b>
60926	<b>Inlet tube with swivel assembly - 1 pack</b>
60927	<b>Inlet tube with swivel assembly - 10 pack</b>
16570	<b>Port cap - 2 pack</b>
60014	<b>Swivel assembly - 10 pack</b>
61905	<b>Mask fitting template</b>



Item #	Description
61013	<b>Mirage Micro™ for Kids nasal mask complete system - small</b> Includes small cushion and small headgear
16333	<b>Mirage Micro™ nasal mask complete system - small</b> Includes small cushion and medium headgear
16334	<b>Mirage Micro nasal mask complete system - medium and large</b> Includes two cushions (medium and large) and medium headgear
16335	<b>Mirage Micro nasal mask complete system - large wide and x-large</b> Includes two cushions (large wide and x-large) and medium headgear
CUSHIONS	
16388	<b>Small cushion</b>
16389	<b>Medium cushion</b>
16390	<b>Large cushion</b>
16391	<b>Large wide cushion</b>
16392	<b>X-large cushion</b>
FRAME	
16378	<b>Mirage Micro and Mirage Micro for Kids nasal mask frame system with small cushion - no headgear</b>
16379	<b>Mirage Micro nasal mask frame system with medium cushion - no headgear</b>
16380	<b>Mirage Micro nasal mask frame system with large cushion - no headgear</b>
16381	<b>Mirage Micro nasal mask frame system with large wide cushion -no headgear</b>
16382	<b>Mirage Micro nasal mask frame system with x-large cushion -no headgear</b>
16383	<b>Mirage Micro and Mirage Micro for Kids nasal mask complete frame assembly, small - no cushion, no headgear</b>
16384	<b>Mirage Micro nasal mask complete frame assembly, standard - no cushion, no headgear</b>
16394	<b>Mask frame only, standard</b>
16395	<b>Mask frame only, small</b>
HEADGEAR	
16117	<b>Universal headgear for Ultra Mirage II, Micro, Activa, Activa LT, SoftGel, and Quattro - medium</b>
16118	<b>Universal headgear for Ultra Mirage II, Micro, Micro for Kids, Activa, Activa LT, SoftGel, Ultra Mirage Full Face Mask and Quattro - small</b>
16119	<b>Universal headgear for Ultra Mirage II, Micro, Activa, Activa LT, SoftGel, Ultra Mirage Full Face Mask and Quattro - large</b>

**Item #** **ACCESSORIES**

16399	<b>Elbow assembly - includes elbow and swivel</b>
60123	<b>Forehead pad - 1 pack</b>
60124	<b>Forehead pad - 10 pack</b>
16393	<b>Forehead support</b>
16569	<b>Headgear clip - 2 pack</b>
16734	<b>Headgear clip - 10 pack</b>
16570	<b>Port cap - 2 pack</b>
16565	<b>Quick release swivel</b>
61289	<b>Forehead Dial</b>
16387	<b>Elbow</b>



Item #	Description
16548	<b>Ultra Mirage™ II nasal mask complete system - standard</b> Includes cushion and headgear
16549	<b>Ultra Mirage II nasal mask complete system - large</b> Includes cushion and headgear
16550	<b>Ultra Mirage II nasal mask complete system - shallow</b> Includes cushion and headgear
16577	<b>Ultra Mirage II nasal mask complete system - shallow wide</b> Includes cushion and headgear
<b>CUSHIONS</b>	
16556	<b>Standard cushion</b>
16557	<b>Large cushion</b>
16558	<b>Shallow cushion</b>
16735	<b>Shallow wide cushion</b>
<b>FRAME</b>	
16726	<b>Ultra Mirage II nasal mask frame system with standard cushion - no headgear</b>
16727	<b>Ultra Mirage II nasal mask frame system with shallow cushion - no headgear</b>
16728	<b>Ultra Mirage II nasal mask frame system with large cushion - no headgear</b>
16736	<b>Ultra Mirage II nasal mask frame system with shallow wide cushion - no headgear</b>
16594	<b>Ultra Mirage II nasal mask complete frame assembly - no cushion, no headgear</b>
<b>HEADGEAR</b>	
16117	<b>Universal headgear for Ultra Mirage II, Micro, Activa, Activa LT, SoftGel, and Quattro - medium</b>
16118	<b>Universal headgear for Ultra Mirage II, Micro, Micro for Kids, Activa, Activa LT, SoftGel, Ultra Mirage Full Face Mask and Quattro - small</b>
16119	<b>Universal headgear for Ultra Mirage II, Micro, Activa, Activa LT, SoftGel, Ultra Mirage Full Face Mask and Quattro - large</b>
<b>ACCESSORIES</b>	
16563	<b>Cushion clip</b>
16566	<b>Elbow assembly - includes elbow, quick release swivel and vent cover</b>
16572	<b>Elbow clip - 1 pack</b>
16573	<b>Elbow clip - 10 pack</b>
60123	<b>Forehead pad - 1 pack</b>
60124	<b>Forehead pad - 10 pack</b>
16590	<b>Forehead support with pads</b>
16569	<b>Headgear clip - 2 pack</b>
16734	<b>Headgear clip - 10 pack</b>
16570	<b>Port cap - 2 pack</b>
16565	<b>Quick release swivel</b>
16575	<b>Vent cover</b>
61905	<b>Mask fitting template</b>





Item #	Description
61030	<b>Pixi™ Pediatric nasal mask complete system</b> Includes cushion and headgear
<b>CUSHIONS</b>	
61033	<b>Pixi Pediatric cushion</b>
<b>FRAME SYSTEM</b>	
61035	<b>Pixi Pediatric frame system - no headgear</b>
<b>HEADGEAR</b>	
61034	<b>Pixi Pediatric headgear</b>



The Mirage Kidsta™ nasal mask is intended for multi-patient use for children aged seven and older or weighing more than 40 lbs, prescribed continuous positive airway pressure (CPAP) or bilevel therapy in hospital, clinic and/or home environments. The Mirage Kidsta nasal mask is intended to be used with the S9 VPAP™ ST-A.

<b>Item #</b>	<b>Description</b>
60000	<b>Mirage Vista™ nasal mask complete system - standard</b> Includes cushion and headgear
60001	<b>Mirage Vista nasal mask complete system - deep</b> Includes cushion and headgear
61011	<b>Mirage Kidsta™ nasal mask complete system - small</b> Includes cushion and headgear
<b>CUSHIONS</b>	
60921	<b>Mirage Vista Standard cushion</b>
60922	<b>Mirage Vista Deep cushion</b>
60013	<b>Mirage Kidsta Small Cushion</b>
<b>FRAME</b>	
60929	<b>Mirage Vista nasal mask frame system with standard cushion - no headgear</b>
60930	<b>Mirage Vista nasal mask frame system with deep cushion - no headgear</b>
60928	<b>Mirage Vista nasal mask complete frame assembly- no cushion, no headgear</b>
60924	<b>Mirage Vista and Mirage Kidsta Mask frame - one size fits all, mask frame only</b>
<b>HEADGEAR</b>	
60918	<b>Mirage Vista Headgear assembly</b>
61020	<b>Mirage Kidsta Headgear assembly</b>
<b>ACCESSORIES</b>	
60919	<b>Mirage Vista Headgear clips - 2 pack</b>
60920	<b>Mirage Vista Headgear clips - 10 pack</b>
61025	<b>Mirage Kidsta Headgear clips - 2 pack</b>
61026	<b>Mirage Kidsta Headgear clips - 10 pack</b>
60926	<b>Inlet tube with swivel assembly - 1 pack</b>
60927	<b>Inlet tube with swivel assembly - 10 pack</b>
16570	<b>Port caps - 2 pack</b>
60014	<b>Swivel assembly - 10 pack</b>
60940	<b>2D fitting template</b>

<b>Item #</b>	<b>Description</b>
16018	Infant Mask System without sensor
<b>CUSHIONS</b>	
16909	Infant cushion Series 5 - large
16910	Infant cushion Series 5 - small
<b>FRAME</b>	
16090	Mask frame
<b>HEADGEAR</b>	
16112	Rescap headgear
<b>ACCESSORIES</b>	
16908	Air tubing Series 5



Item #	Description
61700	<b>Quattro™ FX full face mask complete system - small</b> Includes cushion and medium/standard headgear
61701	<b>Quattro™ FX full face mask complete system - medium</b> Includes cushion and medium/standard headgear
61702	<b>Quattro™ FX full face mask complete system - large</b> Includes cushion and medium/standard headgear
<b>CUSHIONS</b>	
61721	<b>Small cushion</b>
61722	<b>Medium cushion</b>
61723	<b>Large cushion</b>
<b>FRAME SYSTEM</b>	
61724	<b>Quattro FX full face mask frame system with small cushion - no headgear</b>
61725	<b>Quattro FX full face mask frame system with medium cushion - no headgear</b>
61726	<b>Quattro FX full face mask frame system with large cushion - no headgear</b>
<b>HEADGEAR</b>	
61733	<b>Quattro FX full face mask headgear - medium/standard</b>
61734	<b>Quattro FX full face mask headgear - small</b>
61738	<b>Quattro FX full face mask headgear - large</b>



Item #	Description
62501	<b>Quattro™ FX for Her full face mask complete system - small</b> Includes cushion and medium/standard headgear
62502	<b>Quattro™ FX for Her full face mask complete system - medium</b> Includes cushion and medium/standard headgear
<b>CUSHIONS</b>	
61721	<b>Small cushion</b>
61722	<b>Medium cushion</b>
61723	<b>Large cushion</b>
<b>FRAME SYSTEM</b>	
62513	<b>Quattro FX for Her full face mask frame system with small cushion - no headgear</b>
62514	<b>Quattro FX for Her full face mask frame system with medium cushion - no headgear</b>
<b>HEADGEAR</b>	
62506	<b>Quattro FX for Her full face mask headgear (light pink) - small</b>
62507	<b>Quattro FX for Her full face mask headgear (light pink) - standard</b>



Item #	Description
61200	<b>Mirage Quattro™ full face mask complete system - x-small</b> Includes cushion and medium headgear
61201	<b>Mirage Quattro full face mask complete system - small</b> Includes cushion and medium headgear
61202	<b>Mirage Quattro full face mask complete system - medium</b> Includes cushion and medium headgear
61203	<b>Mirage Quattro full face mask complete system - large</b> Includes cushion and medium headgear
<b>CUSHIONS</b>	
61290	<b>X-small cushion with clip</b>
61291	<b>Small cushion with clip</b>
61292	<b>Medium cushion with clip</b>
61293	<b>Large cushion with clip</b>
<b>FRAME</b>	
61260	<b>Mirage Quattro full face mask frame system with x-small cushion - no headgear</b>
61261	<b>Mirage Quattro full face mask frame system with small cushion - no headgear</b>
61262	<b>Mirage Quattro full face mask frame system with medium cushion - no headgear</b>
61263	<b>Mirage Quattro full face mask frame system with large cushion - no headgear</b>
61264	<b>Mirage Quattro full face mask complete frame assembly, x-small - no cushion, no headgear</b>
61265	<b>Mirage Quattro full face mask complete frame assembly, small - no cushion, no headgear</b>
61266	<b>Mirage Quattro full face mask complete frame assembly, medium - no cushion, no headgear</b>
61267	<b>Mirage Quattro full face mask complete frame assembly, large - no cushion, no headgear</b>
61270	<b>Mask frame only, x-small</b>
61271	<b>Mask frame only, small</b>
61272	<b>Mask frame only, medium</b>
61273	<b>Mask frame only, large</b>



<b>HEADGEAR</b>	
16117	Universal headgear for Ultra Mirage II, Micro, Activa, Activa LT, SoftGel, and Quattro - medium
16118	Universal headgear for Ultra Mirage II, Micro, Micro for Kids, Activa, Activa LT, SoftGel, Ultra Mirage Full Face Mask and Quattro - small
16119	Universal headgear for Ultra Mirage II, Micro, Activa, Activa LT, SoftGel, Ultra Mirage Full Face Mask and Quattro - large
<b>ACCESSORIES</b>	
61286	Valve and clip assembly - 1 pack
61287	Valve and clip assembly - 10 pack
61274	Cushion clip, x-small
61275	Cushion clip, small
61276	Cushion clip, medium
61277	Cushion clip, large
61294	Elbow assembly - includes elbow, valve, clip and swivel
61282	Elbow
60123	Forehead pad - 1 pack
60124	Forehead pad - 10 pack
61288	Forehead support
61289	Forehead dial
16565	Quick release swivel
60115	Headgear clips - 2 pack
60116	Headgear clips - 10 pack
61295	Ports cap - 2 pack
61296	Ports cap - 10 pack
61905	Fitting template



Item #	Description
60600	Ultra Mirage™ full face mask complete system - small standard
60601	Ultra Mirage full face mask complete system - small shallow
60602	Ultra Mirage full face mask complete system - medium standard
60603	Ultra Mirage full face mask complete system - medium shallow
60604	Ultra Mirage full face mask complete system - large standard
60605	Ultra Mirage full face mask complete system - large shallow
CUSHIONS	
16604	Small standard cushion
16671	Small shallow cushion
16605	Medium standard cushion
16672	Medium shallow cushion
16606	Large standard cushion
16673	Large shallow cushion
FRAME	
60620	Ultra Mirage full face mask frame system with small standard cushion - no headgear
60621	Ultra Mirage full face mask frame system with small shallow cushion - no headgear
60622	Ultra Mirage full face mask frame system with medium standard cushion - no headgear
60623	Ultra Mirage full face mask frame system with medium shallow cushion - no headgear
60624	Ultra Mirage full face mask frame system with large standard cushion - no headgear
60625	Ultra Mirage full face mask frame system with large shallow cushion - no headgear
60666	Ultra Mirage full face mask complete frame assembly, small - no cushion, no headgear
60667	Ultra Mirage full face mask complete frame assembly, medium - no cushion, no headgear
60668	Ultra Mirage full face mask complete frame assembly, large - no cushion, no headgear
60670	Mask frame only, small
60671	Mask frame only, medium
60672	Mask frame only, large





<b>HEADGEAR</b>	
60674	Headgear - includes headgear clips, one size fits all
16118	Universal headgear for Ultra Mirage II, Micro, Micro for Kids, Activa, Activa LT, SoftGel, Ultra Mirage Full Face Mask and Quattro - small
16119	Universal headgear for Ultra Mirage II, Micro, Activa, Activa LT, SoftGel, Ultra Mirage Full Face Mask and Quattro - large
<b>ACCESSORIES</b>	
16613	Anti-asphyxia valve - includes elbow retainer, valve frame, valve membrane, valve elbow, swivel clip and swivel
16617	Anti-asphyxia valve membrane - 1 pack
16618	Anti-asphyxia valve membrane - 10 pack
16674	Cushion clip, small
16675	Cushion clip, medium
16676	Cushion clip, large
16637	Elbow and swivel clip kit - includes five elbow clips and five quick release clips
60683	Elbow clip
60123	Forehead pad - 1 pack
60124	Forehead pad - 10 pack
60673	Forehead support - 1 pack
60115	Headgear clips - 2 pack
60116	Headgear clips - 10 pack
60681	Luer locks - 2 pack
60682	Luer locks - 10 pack
16689	Swivel only
60680	Swivel clip
16624	Valve elbow
16981	3D fitting template



<b>Item #</b>	<b>Description</b>
61300	<b>Mirage Liberty™ full face mask complete system - small</b> Includes cushion, pillows and headgear
61301	<b>Mirage Liberty full face mask complete system - large</b> Includes cushion, pillows and headgear
<b>CUSHIONS</b>	
61330	<b>Small mouth cushion</b>
61331	<b>Large mouth cushion</b>
<b>PILLOW SLEEVES</b>	
61333	<b>Dual-wall nasal pillows - small</b>
61334	<b>Dual-wall nasal pillows - medium</b>
61335	<b>Dual-wall nasal pillows - large</b>
<b>FRAME</b>	
61318	<b>Mirage Liberty full face mask frame system with small cushion and small pillows - no headgear</b>
61319	<b>Mirage Liberty full face mask frame system with small cushion and medium pillows - no headgear</b>
61320	<b>Mirage Liberty full face mask frame system with small cushion and large pillows - no headgear</b>
61321	<b>Mirage Liberty full face mask frame system with large cushion and small pillows - no headgear</b>
61322	<b>Mirage Liberty full face mask frame system with large cushion and medium pillows - no headgear</b>
61323	<b>Mirage Liberty full face mask frame system with large cushion and large pillows - no headgear</b>
61325	<b>Mirage Liberty full face mask complete frame assembly, small - no cushion, no pillows, no headgear</b>
61326	<b>Mirage Liberty full face mask complete frame assembly, large - no cushion, no pillows, no headgear</b>
61337	<b>Mask frame only, small</b>
61338	<b>Mask frame only, large</b>



<b>HEADGEAR</b>	
61348	Headgear - small
61349	Headgear - standard
<b>ACCESSORIES</b>	
61342	Valve and clip assembly - 1 pack
61343	Valve and clip assembly - 10 pack
61340	Elbow assembly - includes elbow, valve, clip and swivel
61352	Upper clips (left & right)
61353	Lower clips - 2 pack
16570	Ports cap - 2 pack
60926	Inlet tube with swivel assembly - 1 pack
60927	Inlet tube with swivel assembly - 10 pack
60014	Swivel assembly - 10 pack
61905	Fitting template

Item #	Description
16015	<b>Chin Restraint</b> Nasal mask accessory
61914	<b>Gecko™ Nasal Pad: Small - 10 pack (case)</b>
61915	<b>Gecko™ Nasal Pad: Large - 10 pack (case)</b>
16120	<b>Camouflage Headgear</b> Medium - Universal Headgear for Ultra Mirage II, Micro, Activa, Activa LT, SoftGel and Quattro masks
16121	<b>Light Blue Headgear</b> Medium - Universal Headgear for Ultra Mirage II, Micro, Activa, Activa LT, SoftGel and Quattro masks
16122	<b>Pink Headgear</b> Medium - Universal Headgear for Ultra Mirage II, Micro, Activa, Activa LT, SoftGel and Quattro masks
16123	<b>Leopard Headgear</b> Medium - Universal Headgear for Ultra Mirage II, Micro, Activa, Activa LT, SoftGel and Quattro masks



The S9 system sets the new industry standard for comfortable sleep apnea treatment. Equipped with sophisticated new features that showcase ResMed's innovation and dedication to patient compliance, the stylish S9 system delivers everything you and your patients need - and more.

Item #	Description	ADP Billing
36005	<b>S9 AutoSet™</b> Features include the new Enhanced AutoSet algorithm with central sleep apnea detection. Superior APAP comfort with whisper quiet operation and dramatically reduced conducted noise levels. Also includes EPR with Easy-Breathe technology; available in AutoSet mode as well as CPAP mode. Includes usage-based compliance and efficacy downloads for accurate patient monitoring and a standard data card for easy two-way data transfer. Optional accessories support wired and wireless compliance data transfer and integrated pulse oximetry monitoring.	REARJ077D
36015	<b>S9 AutoSet &amp; H5i™</b> S9 AutoSet packaged with the H5i humidifier (Item #36900).	REARJ077D
36075	<b>S9 AutoSet &amp; H5i for Her (light pink)</b> S9 AutoSet for Her packaged with the H5i for Her humidifier.	REARJ077D
36025	<b>S9 AutoSet &amp; Climate Control Kit</b> S9 AutoSet with the H5i humidifier and the ClimateLine™ Tubing (Item #36995)	RENLV045M
36003	<b>S9 Elite™</b> Features include new central sleep apnea detection, whisper quiet operation and dramatically reduced conducted noise levels. Also includes EPR with Easy-Breathe technology; Includes usage-based compliance and efficacy downloads for accurate patient monitoring and a standard data card for easy two-way data transfer. Optional accessories support wired and wireless compliance data transfer and integrated pulse oximetry monitoring.	RENLV045M
36013	<b>S9 Elite &amp; H5i</b> S9 Elite packaged with the H5i humidifier (Item #36900).	RENLV045M
36023	<b>S9 Elite &amp; Climate Control Kit</b> S9 Elite with the H5i humidifier and the ClimateLine Tubing (Item #36995).	RENLV045M



The S9 system sets the new industry standard for comfortable sleep apnea treatment. Equipped with sophisticated new features that showcase ResMed's innovation and dedication to patient compliance, the stylish S9 system delivers everything you and your patients need - and more.

Item #	Description	ADP Billing
36002	<b>S9 Escape Auto™</b> Features include superior APAP mode based upon the gold-standard AutoSet Algorithm providing comfort with whisper quiet operation and dramatically reduced conducted noise levels. Also includes EPR with Easy-Breathe technology; available in AutoSet mode as well as CPAP mode. Includes usage-based compliance and AHI data for accurate patient monitoring, along with a standard data card for easy two-way data transfer. Optional accessories support wired and wireless compliance data transfer.	
36012	<b>S9 Escape Auto &amp; H5i™</b> S9 Escape Auto packaged with the H5i humidifier (Item #36900).	
36022	<b>S9 Escape Auto &amp; Climate Control Kit</b> S9 Escape Auto with the H5i humidifier and the ClimateLine Tubing (Item #36995)	
36001	<b>S9 Escape™</b> Features include whisper quiet operation and dramatically reduced conducted noise levels. Also includes EPR with Easy-Breathe technology; Includes usage-based compliance downloads for accurate patient monitoring along with a standard data card for easy two-way data transfer. Optional accessories support wired and wireless compliance data transfer.	RENLV045N
36011	<b>S9 Escape &amp; H5i</b> S9 Escape packaged with the H5i humidifier (Item #36900).	RENLV045N
36021	<b>S9 Escape &amp; Climate Control Kit</b> S9 Escape with the H5i humidifier and the ClimateLine Tubing (Item #36995)	RENLV045N

Item #	Description	ADP Billing
36004	<b>S9 VPAP™ S</b> A bilevel device featuring CPAP and S therapy modes. The S9 VPAP S includes the comfort option of Easy-Breathe waveform and central sleep apnea (CSA) detection in S mode. Optional oximetry monitoring and wireless data compliance transfer.	PS0000106
36014	<b>S9 VPAP S &amp; H5i</b> S9 VPAP S packaged with the H5i humidifier.	PS0000106
36024	<b>S9 VPAP S &amp; Climate Control Kit</b> S9 VPAP S with the H5i humidifier and the ClimateLine tubing.	PS0000106
36006	<b>S9 VPAP™ Auto</b> An auto-adjusting bilevel device featuring CPAP, S and VAuto therapy modes. The VAuto, auto-adjusting bilevel mode uses ResMed's enhanced AutoSet algorithm to preemptively respond to events and detect central sleep apnea (CSA). Easy-Breathe waveform is standard in VAuto mode and a comfort option in S mode. Optional oximetry monitoring and wireless compliance data transfer.	PS0000107
36016	<b>S9 VPAP Auto &amp; H5i</b> S9 VPAP Auto packaged with the H5i humidifier.	PS0000107
36026	<b>S9 VPAP Auto &amp; Climate Control Kit</b> S9 VPAP Auto with the H5i humidifier and the ClimateLine Tubing.	PS0000107
36008	<b>S9 VPAP™ ST</b> A fully-integrated bilevel device that delivers comfortable and effective therapy with exceptional patient-ventilator synchrony to patients requiring backup rate. It also features optional integrated oximetry and a sleep quality indicator, which allows patients to view their therapy progress on a daily basis.	
36018	<b>S9 VPAP ST &amp; H5i</b> S9 VPAP ST packaged with the H5i humidifier.	
36028	<b>S9 VPAP ST &amp; Climate Control Kit</b> S9 VPAP ST with the H5i humidifier and the ClimateLine tubing.	
36037	<b>S9 VPAP™ Adapt</b> An adaptive servo-ventilator specifically designed to treat central sleep apnea (CSA) in all its forms, including mixed events now with optional ASV Auto mode that stabilizes the upper airway to treat obstructive apneas while maintaining stable minute ventilation and treating all forms of periodic breathing. It also features optional integrated oximetry and a sleep quality indicator, which allows patients to view their therapy progress on a daily basis.	
36047	<b>S9 VPAP Adapt &amp; H5i™</b> S9 VPAP Adapt packaged with the H5i humidifier.	
36057	<b>S9 VPAP Adapt &amp; Climate Control Kit</b> S9 VPAP Adapt with the H5i humidifier and the ClimateLine tubing.	

<b>HUMIDIFIERS</b>	
36900	H5i - for use with S9 Flow Generators
36901	H5i with cleanable water tub - for use with S9 Flow Generators
<b>H5i HUMIDIFIER ACCESSORIES</b>	
36995	H5i ClimateLine Tubing - 1 pack
36997	H5i ClimateLine MAX™ Tubing - 1 pack
36803	H5i Standard Water Tub - is an openable tub for easy cleaning - replaced Item #36802 H5i Water Tub (sealed)
36800	H5i Cleanable Water Tub
36891	H5i Flip Lid
36892	H5i Flip Lid Seal
36894	H5i Air Outlet and Outlet Clip
<b>S9 FLOW GENERATOR ACCESSORIES</b>	
36810	SlimLine™ Tubing - for use with S9 and Stellar
36811	SlimLine Tubing Wrap
24947	ResMed Power Station II Battery - requires S9 90W Power Supply to charge
36821	S9 90W Power Supply Unit Plus Power Cord - can be used with or without humidifier - for use with S9 and ResMed Power Station II Battery
36920	S9 30W Power Supply Unit - for use with device only, will not support humidifier
36860	S9 Travel Bag
36859	S9 Filter Cover
36970	DC/DC converter - 24V/90W output
<b>S9 FILTERS</b>	
36850	S9 Standard Filter - 1 pack
36851	S9 Standard Filter - 2 pack
36852	S9 Standard Filter - 12 pack
36853	S9 Standard Filter - 50 pack
36855	S9 Hypo-allergenic Filter - 1 pack
36856	S9 Hypo-allergenic Filter - 2 pack
36857	S9 Hypo-allergenic Filter - 12 pack
36858	S9 Hypo-allergenic Filter - 50 pack





Stellar Ventilators	
24157	<b>Stellar™ 100 Adult and Pediatric Ventilator</b> Ultra-premium non-invasive ventilator with invasive option for adult and pediatric (13 kg) ventilation. Pressure from 3-40 cm H <sub>2</sub> O, PS to 37 cm H <sub>2</sub> O, and internal battery. FiO <sub>2</sub> and SpO <sub>2</sub> monitoring (optional)
24140	<b>Stellar™ 150 Adult and Pediatric Ventilator</b> Ultra-premium non-invasive ventilator with invasive option for adult and pediatric (13 kg) ventilation. Comprehensive ventilation modes including iVAPS (Intelligent Volume Assured Pressure Support). Pressure from 3-40 cm H <sub>2</sub> O, PS to 37 cm H <sub>2</sub> O, and internal battery. FiO <sub>2</sub> and SpO <sub>2</sub> monitoring (optional).
Stellar Humidifier and Accessories	
24927	<b>H4i Heated Humidifier for Stellar</b>
36810	<b>SlimLine tubing - for use with S9 and Stellar</b>
24977	<b>Stellar Mobility bag</b>
22305	<b>XPOD Oximeter</b>
1431002	<b>Oximeter Flex Sensor (replaces obsolete item #16936) - includes sensor and item #70276 (single use adhesive strips, 25 pack)</b>
24955	<b>FiO<sub>2</sub> Monitoring Kit</b>
CAP008617	<b>FiO<sub>2</sub> Monitoring Sensor</b>
24976	<b>ResMed Leak Port</b>
24933	<b>Stellar Filter - hypoallergen - 1 pack</b>
24934	<b>Stellar Filter - hypoallergen - 12 pack</b>
24905	<b>Standard stand (without top plate) - must purchase top plate separately to work with PAP devices (see Item #24913)</b>
24913	<b>VPAP Series top plate - for use with S7, S8, S9 and Stellar</b>
249007	<b>Hospital Stand with Stellar Top Plate</b>
24917	<b>Premium circuit support arm</b>
24919	<b>Tube clip - double</b>
24920	<b>Tube clip - single</b>
24966	<b>Antibacterial Filter - 50 pack</b>
24967	<b>HME filter - Adult 50 pack</b>
24968	<b>HME filter - Pediatric Small Angled 50 pack</b>
24969	<b>HME filter - Pediatric 50 pack</b>
14948	<b>Grey Tubing - cuffed</b>
24972	<b>AC Power Cord Figure 8</b>
24935	<b>ResMed USB Data Storage Stick</b>
24971	<b>O<sub>2</sub> Coupling Insert</b>



Item #	S9 VPAP ST-A and S9 VPAP ST
36039	<p><b>S9 VPAP™ ST-A</b> A fully-integrated noninvasive ventilation device designed to treat non-dependent patients requiring treatment pressure up to 30 cm H<sub>2</sub>O with the security of alarms and a respiratory backup rate up to 50 BPM. It includes iVAPS mode, treating patients with changing respiratory needs. It also features optional integrated oximetry and a sleep quality indicator, which allows patients to view their therapy progress on a daily basis.</p>
36049	<p><b>S9 VPAP ST-A &amp; H5i</b> S9 VPAP ST-A packaged with the H5i humidifier (Item #36900).</p>
36059	<p><b>S9 VPAP ST-A &amp; Climate Control Kit</b> S9 VPAP ST-A with the H5i humidifier and the ClimateLine MAX™ tubing (Item #36997).</p>
36008	<p><b>S9 VPAP™ ST</b> A fully-integrated bilevel device that delivers comfortable and effective therapy with exceptional patient-ventilator synchrony to patients requiring backup rate. It also features optional integrated oximetry and a sleep quality indicator, which allows patients to view their therapy progress on a daily basis.</p>
36018	<p><b>S9 VPAP ST &amp; H5i</b> S9 VPAP ST packaged with the H5i humidifier (Item #36900).</p>
36028	<p><b>S9 VPAP ST &amp; Climate Control Kit</b> S9 VPAP ST with the H5i humidifier and the ClimateLine tubing (Item #36995).</p>



<b>COMPLIANCE SOFTWARE</b>	
31315	<b>ResScan software</b> Clinical management software for use with all ResMed flow generator products.
N/A	<b>EasyCare Online</b> Compliance Management Solution - please contact ResMed Customer Service for account setup.
<b>S9 DATA ACCESSORIES</b>	
36840	<b>Protective Folder with SD card - 10 pack</b>
36845	<b>Protective Folder without SD card - 10 pack</b>
36931	<b>SD Card Reader</b>
36950	<b>S9 USB Adapter</b>
36930	<b>S9 Wireless Module: Transmits summary data wirelessly for use with EasyCare Online</b>
<b>S8 DATA ACCESSORIES</b>	
22903	<b>ResScan data cards with protective folder - 10 pack</b>
22906	<b>Protective Folder (no data card) - 10 pack</b>
22207	<b>ResScan card reader</b>
22212	<b>SmartMedia card reader II - USB</b>
30942	<b>SmartMedia card - 2 pack</b>
22201	<b>Serial port adapter</b>
22203	<b>ResScan USB serial adapter</b>
22915	<b>Data module cover</b>
17952	<b>Null modem cable: for use with VPAP III, AutoSet Spirit™, AutoSet Respond™, S7™ Elite CPAP, S8 CPAP and AutoSet T flow generators</b>
<b>S9 OXIMETRY</b>	
36940	<b>S9 Oximeter Adapter</b>
369100	<b>S9 Complete Oximetry Kit</b> (includes Item #36940 and Item #36942)
<b>S8 OXIMETRY</b>	
22206	<b>S8 ResLink™</b> Includes 32 MB SmartMedia card, user's manual and mailing bag.
222102	<b>S8 Complete Oximetry Kit</b> (includes Item #22206 and Item #1430050)
<b>S7 OXIMETRY</b>	
30924	<b>S7 ResLink PMM: patient management module.</b> Includes 32 MB SmartMedia card, user's manual and mailing bag.



Item #	OXIMETRY ACCESSORIES
36942	<b>S9 Oximetry Accessory Kit</b> (includes reusable soft sensor, 3M cable, sensor clip, XPOD)
707560	<b>Reusable soft sensor (medium) - 3M cable</b>
1431002	<b>Disposable flex sensor - 3M cable</b> includes sensor and item #70276 (adhesive strips, single use - 25 pk)
22961	<b>ResLink converter cable</b> This converter cable is for use with ApneaLink™ XPOD to connect the oximeter port for either the S7 ResLink or S8 ResLink modules.
36941	<b>S9 Oximeter converter cable</b> This converter cable is for use with existing ResLink XPOD to connect with the S9 Oximeter module.



Item #	Description
22328	<p><b>ApneaLink™ Plus Complete Kit</b> A portable home sleep test diagnostic device that records up to four channels of information: respiratory effort, pulse, oxygen saturation and nasal flow. Kit includes ApneaLink Plus recorder, software CD, quick software setup guide, carry bag, EasySense respiratory effort sensor, 2 reusable belts, USB download cable, 2 AA batteries, 3 nasal cannulas, XPOD oximeter, 3 single use oximetry sensors and sensor clip.</p>
22319	<p><b>ApneaLink Plus Basic Kit</b> A portable home sleep test diagnostic device that records up to four channels of information: respiratory effort, pulse, oxygen saturation and nasal flow. Kit includes ApneaLink Plus recorder, software CD, quick software setup guide, carry bag, EasySense respiratory effort sensor, 2 reusable belts, USB download cable, 2 AA batteries and 1 nasal cannula.</p>
22302	<p><b>ApneaLink</b> A sleep screening device with up to three channels (with oximetry) used to objectively assist in identifying the risk for sleep-disordered breathing. Includes ApneaLink software CD, operating instructions, carry bag, 3 nasal cannulas, single use belt, USB download cable and 2 AA batteries.</p>
22304	<p><b>ApneaLink Oximetry Kit</b> Provides oximetry and pulse information to the ApneaLink and ApneaLink Plus recorder devices. Uses single-use or reusable probes. Must be used in conjunction with ApneaLink and ApneaLink Plus only. Kit includes oximeter (XPOD), oximeter sensor clip and 3 single use sensors.</p>



<b>ACCESSORIES AND CONSUMMABLES</b>	
22339	<b>ApneaLink Plus and ApneaLink Software CD 9.3 (compatible with all devices)</b>
22333	<b>ApneaLink Plus EasySense Respiratory Effort Sensor (for use with ApneaLink Plus ONLY)</b>
<b>NASAL CANNULAS</b>	
70388	<b>ApneaLink standard cannula - 25 pack</b>
<b>BELTS</b>	
629052	<b>Reusable belt - for use with ApneaLink and ApneaLink Plus</b>
70406	<b>Single use belt - 10 pack (for use with ApneaLink ONLY)</b>
<b>MISCELLANEOUS</b>	
629051	<b>USB download cable</b>
22309	<b>Rechargeable battery kit - includes battery charger and 4 AA rechargeable batteries</b>
<b>OXIMETRY COMPONENTS</b>	
22306	<b>Oximeter sensor clip</b>
70413	<b>Reusable soft sensor - large</b>
70567	<b>Reusable soft sensor - small</b>
70726	<b>Reusable flex sensor kit, 36" includes sensor and Item #70276</b>
70276	<b>Adhesive strips, single use - for use with Item #1431002 or Item #70726 - 25 pack</b>
22338	<b>Single use sensor, adult, 36" (1 meter) (ResMed Standard) - 24 pack. Replaces Item #22316</b>
22337	<b>Single use sensor, 36" (1 meter) (value line) - 24 pack. Replaces Item #70412</b>
36955	<b>S9 Embletta Adapter: to be used with Embletta HST device. Requires separate purchase of 3-in-1 cable from Embla.</b>

<b>POWER SUPPLIES</b>	
70199	<b>Battery adapter, 12 volt</b>
26932	<b>DC converter - for VPAP ST, VPAP Auto 25 and VPAP S</b>
22015	<b>DC converter - for VPAP III ST-A</b>
33942	<b>DC-12 converter - for S8 CPAP unit and VPAP Auto</b>
36830	<b>S9 Power cord</b> This black power cord can be used for S9, S8 and C-Series Tango Units. This Item does not include the main power supply. (replaced Item #33938)
26905	<b>Power cord - for VPAP Adapt SV and VPAP III ST-A</b>
33948	<b>S8 connector plug - 10 pack</b>
<b>RESMED POWER STATION II</b>	
24947	<b>ResMed Power Station II - Battery pack (requires S9 90W Power Supply to charge)</b>
24923	<b>ResMed Power Station II - Battery pack, PSU adaptor and bag (requires S9 90W power supply to charge)</b>
24959	<b>DC cable, Stellar and VPAP III</b>
24981	<b>DC cable, Stellar, 0.5m</b>
24961	<b>DC cable, S9 Series</b>
24960	<b>DC cable, S8 VPAP</b>
24946	<b>PSU adaptor, ResMed Power Station II</b>
36821	<b>90W Power Supply Unit Plus Power Cord - can be used with or without humidifier - for use with S9 and ResMed Power Station II Battery</b>
24964	<b>Battery coupler kit, ResMed Power Station II</b>
24965	<b>Battery coupler, ResMed Power Station II</b>
24952	<b>Carry bag, ResMed Power Station II</b>
24943	<b>Velcro strap, ResMed Power Station II</b>
249002	<b>ResMed Power Station II kit for Stellar/S7</b>
249003	<b>ResMed Power Station II kit for S8 VPAPs</b>
249004	<b>ResMed Power Station II kit for S9</b>
<b>TUBINGS AND AIR DELIVERY SYSTEMS</b>	
26909	<b>Air delivery system - for AutoSet CS 2, VPAP Adapt SV Enhanced Includes sensor tubing, tube clips, proximal cuff and standard 2 meter air hose</b>
14911	<b>Air tubing connector - 5 pack</b>
14994	<b>Air tubing, standard - light gray, ribbed (6')</b>
14999	<b>Air tubing, clear gray ribbed (9', 2.75 meter) - for use with all S7, S8 products and S9 VPAP only</b>
14986	<b>Air tubing, autoclave - dark gray (6'6", 2 meter)</b>
14922	<b>Air tubing, gray (9'9", 3 meter)</b>
14987	<b>Air tubing, standard - light gray, cuffed (6'6", 2 meter)</b>
26960	<b>Oxygen adapter port - for VPAP ST, VPAP Auto 25 and VPAP S</b>
14915	<b>Oxygen sideport connector</b>
26915	<b>Proximal cuff - for AutoSet CS 2 and VPAP Adapt SV Enhanced</b>
26910	<b>Sensor tubing with Luer lock - for AutoSet CS 2 and VPAP Adapt SV Enhanced</b>
26911	<b>Tube clips - 9 pack - for AutoSet CS 2 and VPAP Adapt SV Enhanced</b>

CARRY BAGS / TRAVEL BAGS / STAND	
33972	<b>S8 Series II premium travel bag (black)</b>
30964	<b>Hospital cart/stand (S7, VPAP III) - top plate included</b>
AIR FILTER COVERS	
30909	<b>S7 Series and VPAP III air filter cover</b>
33923	<b>Filter cover for S8</b>
33959	<b>Blue filter cover for S8</b>
MANOMETERS AND PRESSURE SENSOR ADAPTERS	
30901	<b>Pressure sensor adapter</b> Clear plastic adapter for external humidifier connection with the AutoSet Spirit, AutoSet Respond and S7 Elite
30916	<b>Pressure sensor adapter</b> Rubber adapter for AutoSet Spirit, AutoSet Respond and S7 Elite CPAP



Item #	Description
33915	S8, S8 II and C-Series Tango filters - 1 pack
33973	S8, S8 II and C-Series Tango filters - 2 pack
33974	S8, S8 II and C-Series Tango filters - 12 pack
33918	S8, S8 II and C-Series Tango filters - 50 pack
33919	S8, S8 II and C-Series Tango hypo filters - 1 pack
33975	S8, S8 II and C-Series Tango hypo filters - 2 pack
33976	S8, S8 II and C-Series Tango hypo filters - 12 pack
33922	S8, S8 II and C-Series Tango hypo filters - 50 pack
30919	S7 Series and VPAP III filters - 1 pack
30907	S7 Series and VPAP III filters - 3 pack
30908	S7 Series and VPAP III filters - 50 pack





Item #	Description
26940	<b>H4i - for use with all S8, S8 II and VPAP bilevel products.</b>
26950	<b>H4i with optional cleanable water chamber</b> for use with all S8, S8 II and VPAP bilevel products.
30902	<b>HumidAire 2i™</b> for use with S7™ Lightweight CPAP, S7 Elite CPAP, AutoSet Spirit, AutoSet Respond, VPAP III and VPAP Adapt SV

	Description
26952	<b>H4i water chamber (standard)</b>
26958	<b>H4i water chamber (cleanable)</b>
26953	<b>H3i top cover and seal for H4i water chamber</b>
26954	<b>H4i tub plate seal</b>
26955	<b>H4i cleaning plug</b>
26963	<b>H4i spill guard</b>
33926	<b>HumidAire 3i top cover and seal</b>
33927	<b>HumidAire 3i top cover seal</b>
33947	<b>HumidAire 3i cleaning plug pack</b>
30951	<b>HumidAire 2i tub assembly</b>
30915	<b>HumidAire 2i water chamber o-ring</b>
30914	<b>HumidAire 2i water chamber tub plate</b>
30938	<b>HumidAire 2i water chamber cross member, PC</b>
30961	<b>HumidAire 2i base PC-side catch</b>
30939	<b>HumidAire 2i latch, PC</b>
30955	<b>HumidAire 2i side catches, PC</b>
30956	<b>HumidAire 2iC tub assembly</b>
30912	<b>HumidAire 2i and 2iC water chamber seals</b>
34901	<b>C-Series water chamber</b>
34903	<b>C-Series adapter seal - 10 pack</b>
18920	<b>HumidAire water chamber</b>
18913	<b>HumidAire water chamber base</b>
18911	<b>HumidAire water chamber lid</b>
18912	<b>HumidAire water chamber seal</b>
33963	<b>Tubing wrap (6'6") - reduces rainout in tubing caused by heated humidification</b>

This Home Care Provider Price List has been prepared solely for the use of ResMed's customers. ResMed reserves the right to modify any of the information contained herein at any time without notice, including pricing, product descriptions and configurations, and warranty terms.

All information contained in this Price List is confidential. You agree not to disclose such information to any third party without the prior written approval of ResMed. While ResMed has taken reasonable precautions to ensure the accuracy of the pricing and information presented in this document as of the date of publication, ResMed strongly recommends that you contact customer service and/or a sales representative to confirm the accuracy of the information or to see if other discounts, promotions, or specials may apply.

All purchases of ResMed products by Customer will be governed by ResMed's Standard Terms and Conditions of Sales located on the back of each invoice, unless otherwise agreed upon in separate writing by ResMed.

Please contact ResMed's Customer Service Department at (800) 424-0737 if you have any questions or need information about ResMed's standard terms and conditions.



ResMed Corp San Diego, CA, USA +1 858 836 5000 or 1 800 424 0737 (toll free), ResMed Ltd Bella Vista, NSW, Australia +61 (2) 8884 1000 or 1 800 658 189 (toll free).

Offices in Austria, Brazil, Finland, France, Germany, Hong Kong, Japan, Malaysia, Mexico, Netherlands, New Zealand, Norway, Singapore, Spain, Sweden, Switzerland, United Kingdom (see website [www.resmed.com](http://www.resmed.com) for details). Protected

Activa, AutoScan, AutoSet, AutoSet Respond, AutoSet Spirit, ClimateLine, EPR, Gecko, H4i, HumidAire, Kidsta, Liberty, Mirage, Mirage Activa, Mirage FX, Mirage Kidsta, Mirage Liberty, Mirage Micro, Mirage Quattro, Mirage Swift, Mirage Vista, Pixi, Quattro, ResControl, ResLink, ResScan, ResTraxx, SmartStart, S7 Elite, S7 Lightweight, S8, S8 Compact, S8 Elite, S8 Escape, S9, S9 AutoSet, S9 Elite, S9 Escape, S9 Escape Auto, SlimLine, Stellar, Sullivan, Swift, Tango, TiControl, Ultra Mirage, Vista, VPAP, VPAP Adapt SV, VPAP Auto, VPAP Malibu and VSet are trademarks of ResMed Ltd.

Activa, ApneaLink, AutoSet, Escape, HumidAire, Kidsta, Mirage, Mirage Quattro, SmartStart, Sullivan, Vista and VPAP are registered in U.S. Patent and Trademark Office.

ApneaLink is a trademark of MAP Medizin Technologie GmbH and Breath-O-Prene™ is a trademark of AccuMed. ©2007 ResMed Ltd